

Investigation of extraneous peak

1 Investigation of extraneous peak reveals ethylene oxide in food supplements

2 Abdulla Siddiqui, Ahmed Jaber Al Shareef, & Premanandh Jagadeesan*

3 *Technical Affairs Division, Central Testing Laboratories,*

4 *Quality and Conformity Council*

5 *P.O Box 853, Abu Dhabi, United Arab Emirates*


6

7

8 *Corresponding Author: jpanandh@yahoo.com

9

10 **Graphical abstract**



11

12 The global food supplement market is growing rapidly and so as the concerns on their
13 quality, safety, and product efficacy. Product monitoring is challenging due to lenient
14 regulation across the globe when compared to pharmaceutical counterparts. As part of
15 quality and safety compliance of food supplements, products are being tested randomly in
16 the Emirate of Abu Dhabi, UAE. This study is the part of routine monitoring of
17 supplements during which an extraneous peak was observed and characterized. Analysis
18 was performed using a Gas chromatography coupled to mass spectrometry and headspace
19 extraction sampler. The extraneous peak with signal to noise ratio (S/N) of 140 was
20 identified as ethylene oxide. Since traces were observed in active product formulation, it
21 is worthwhile to monitor the levels in a higher number of supplements to elucidate its
22 significance.

23 Key words: supplement; ethylene oxide; GC-MS-MS; contaminants; extraneous peak;
24 maximum residue limit

25

26 **Introduction**

27 Food is the basic sustenance for life and food supplements are substances
28 intended to complement dietary requirements having nutritional or physiological
29 effects. In general, food supplements comprise of nutritive products such as
30 vitamins, minerals, amino acids, fatty acids, and fibers that are essential for the
31 maintenance of good health. The Food, Drug, and Cosmetic Act of Food and
32 drug administration (FDA) defines a dietary supplement as a product, taken
33 orally, containing a dietary ingredient intended to supplement the diet [1].

34 The global food supplement market has grown over 140 billion USD in recent
35 times and is estimated to have a compound annual growth rate (CAGR) of
36 around 9.0% by 2028 [2]. Thus, increasing use of dietary supplements claiming
37 health benefits raises concerns on their quality and safety in addition to product
38 efficacy. Also, product monitoring is challenging due to lenient regulation across
39 the globe when compared to pharmaceutical counterparts [3].

40 In the United Arab Emirates (UAE), ethnically diverse group of people hailing
41 from different cultural and traditional background live and use of different forms of
42 dietary supplement is not uncommon [4]. According to one study, more than half
43 of the studied population had taken dietary supplements at some point [2].

44 From the regulatory perspective, dietary supplements having classification as
45 “foods” by the Ministry of Health is controlled by food control authority while the
46 borderline products are submitted to the Ministry of Health for classification [5].

47 Food supplements with specific intention to enhance nutritional elements such as

Investigation of extraneous peak

48 vitamins, minerals, fatty acids, amino acids, enzymes, prebiotics, probiotics,
49 collagen, dietary fibers, melatonin, propolis, pollen, herbs or food herbal extracts
50 are covered by UAE.S GSO 2571:2021. As part of quality and safety compliance
51 of food supplements, products are being tested randomly to elucidate their
52 status. This study is the part of routine monitoring of supplements during which
53 an extraneous peak was observed and characterized.

54 **Materials and Methods:**

55 *Sample preparation*

56 The samples were homogenized using a mechanical shaker for 10 minutes and 1
57 gram of the homogenate was added to 50 ml of diluent (acetonitrile) and
58 vortexed for 5 minutes. The mixtures were centrifuged at 5000 rpm for 5 minutes.
59 An aliquot of 5 ml from the centrifuged supernatant was transferred to HS vial
60 and analyzed using GC-MS/MS equipped with dynamic head space sampler.

61 *Standard preparation*

62 Reference standard (RESTEK # 30548) was used to create calibration curves.
63 The final concentration of ethylene oxide matrix matched standards solutions
64 was prepared at 5, 10, 20, 50, 100 and 200 ng/ml levels in acetonitrile.

65 *Analytical condition:*

66 All analyses were performed using a GC-2010 PLUS system coupled to mass
67 spectrometry and headspace extraction sampler "HS" (Auto sampler AOC-6000
68 Plus-Shimadzu, Japan). The residual analytes were separated on a Rxi-624 Sil
69 MS capillary column. Helium gas was used as the carrier having the initial flow
70 rate of 1.5 mL/min followed by a constant linear velocity of 44 cm/s. Since

Investigation of extraneous peak


71 splitless injection is ideal for trace analysis for very low analyte concentrations, it
72 was used for 0.5 min followed by split mode (5:1) The oven temperature gradient
73 program was as follows: hold an initial temperature of 35 °C for 2 min, ramp to
74 235 °C at 20 °C/min, and finally hold for 2 min. The column was reconditioned at
75 200 °C for 2 min to eliminate all impurities co-extracted from the matrix and
76 minimize carry-over effects. The temperatures of the injection port and MS
77 transfer line were 200 and 230 °C, respectively. The GC instrument was directly
78 interfaced to a GCMSTQ8040 triple quadrupole mass spectrometer (Shimadzu).
79 The MS analyses were conducted in positive electron ionization (EI) mode with a
80 filament current at 60 µA. An ionization energy at 70 eV was applied. The ion
81 source temperature was set at 230 °C. For ion fragmentation, argon (grade 5.0)
82 was used as the collision-induced dissociation (CID) gas. Quantifications were
83 performed using multiple reaction monitoring (MRM) transitions. The specific
84 transitions and the optimum collision energies (CEs) for all compounds are listed
85 in Table 1. The data handling and system operations was carried out by GCMS
86 Solution software (version 4.52; Shimadzu).

87 Results and Discussion:

88 In general, nutritional supplement product ingredients are very diverse in their
89 chemical composition and often do not match label claims [6]. Accidental
90 contamination and economically motivated substances in the supplements raise
91 concerns on the product quality especially with herbal and botanical ingredients,
92 [7].

Investigation of extraneous peak

93 The current investigation of dietary supplement was part of routine monitoring for
94 both chemical and microbiological parameters as per UAE.S GSO 2571:2021. In
95 addition, a forensic screening for volatile organics was also employed with mass
96 ranging from 40-500.



97

98 **Fig.1.calibration curve of ethylene oxide standards and chromatogram exhibiting the retention time.**

99

100 The forensic screening revealed an extraneous peak with signal to noise ratio
101 (S/N) of 140 which does not match with any of the solvents used. Upon further
102 scrutiny, this peak was identified as ethylene oxide based on the mass number
103 obtained by GC-MS with EI source (Fig. 1). Analyte identification criteria was
104 based on retention time (RT), and the relative MRM transition intensities which

Investigation of extraneous peak

105 should not vary by more than $\pm 1\%$ and $\pm 20\%$ respectively, relative to a spiked
106 control sample [8]. It is a disinfectant and effective fumigant used to protect
107 against microorganisms and hence used in the sterilization of tools in medicine or
108 in food industries to extend the nutritional value and shelf life. [9].

109 Of the 4 tested samples, 2 of them were found to contain 23 and 44 ppb of
110 ethylene oxide. The maximum residue limit (MRL) for food supplements is set at
111 0.1 mg/kg (100ppb) including the capsules used in supplement formulations. [10].

112 Though the levels are within the prescribed limits and the traces could attribute to
113 the sterilization of gelatin capsule, we used only the active medicinal compound
114 from the capsule and hence it is worthwhile to monitor the levels in a higher
115 number of supplements to elucidate its significance.

116 In conclusion, the present study has identified the extraneous peak in nutritional
117 supplement products as ethylene oxide demonstrating the need to extend the
118 study to a higher number of samples.

119

120 **Acknowledgments**

121 The authors gratefully acknowledge the support and encouragement of the Abu
122 Dhabi Quality and Conformity Council. Gratitude is also extended to the
123 management of Central Testing Laboratories (CTL) for invaluable support.

124

125 **References**

126 1. Dietary Supplement Health and Education Act of 1994. Public Law 103–
127 417, 103rd Congress. <http://www.fda.gov/opacom/laws/dshea.html>.

- 128 2. Jairoun, A. A., Saleh Al-Hemyari, S., Shahwan, M., Saeed, B. Q.,
129 Godman, B., Kurdi, A., Zyoud, S. H., & Ibrahim, B. (2023). Risk
130 Awareness and Demographic Characteristics Associated with the Use of
131 Sexual Enhancement Supplements Among University Staff: A Cross-
132 Sectional Study in the United Arab Emirates. *American journal of men's*
133 *health*, 17(2), 15579883231158531.
- 134 3. Djaoudene O, Romano A, Bradai YD, Zebiri F, Ouchene A, Yousfi Y,
135 Amrane-Abider M, Sahraoui-Remini Y, Madani K. A Global Overview of
136 Dietary Supplements: Regulation, Market Trends, Usage during the
137 COVID-19 Pandemic, and Health Effects. *Nutrients*. 26;15(15):3320.
- 138 4. Radwan, H., Hasan, H., Jaafar, Z., Abbas, N., Rashed Saif, E., Al Kitbi,
139 M., et al. (2022). Diets and dietary supplements used during the COVID-
140 19 pandemic in the United Arab Emirates: a cross-sectional survey. *Saudi*
141 *Pharm. J.* 30, 421–432. doi: 10.1016/j.jsps.2022.01.019.
- 142 5. Code of Practice, 27. (2019). Health, Nutrition and Related Claims. Abu
143 Dhabi Agricultural and food safety Authority. Available at
144 [https://www.adafsa.gov.ae/English/PolicyAndLegislations/Guidelines/Docu](https://www.adafsa.gov.ae/English/PolicyAndLegislations/Guidelines/Documents/Health,%20Nutrition%20and%20Related%20Claims%20Cop.pdf)
145 [ments/Health,%20Nutrition%20and%20Related%20Claims%20Cop.pdf](https://www.adafsa.gov.ae/English/PolicyAndLegislations/Guidelines/Documents/Health,%20Nutrition%20and%20Related%20Claims%20Cop.pdf)
- 146 6. Dwyer, J. T, Bailen R. A, Saldanha L. G, Gahche J. J, Costello R. B, Betz
147 J. M, Davis C. D, Bailey R. L, Potischman N, Ershow A. G, Sorkin B. C,
148 Kuszak A. J, Rio-Avila L, Chang F, Goshorn J, Andrews K. W, Pehrsson
149 P. R, Gusev P. A, Harnly J. M Hardy C. J, Emenaker N. J, Herrick K

Investigation of extraneous peak

- 150 (2018). The Dietary Supplement Label Database: Recent Developments
151 and Applications. *J Nutr*, 148: 1428S-1435.
- 152 7. Martínez-Sanz, J. M., Sospedra, I., Ortiz, C. M., Baladía, E., Gil-Izquierdo,
153 A., & Ortiz-Moncada, R. (2017). Intended or Unintended Doping? A
154 Review of the Presence of Doping Substances in Dietary Supplements
155 Used in Sports. *Nutrients*, 9(10), 1093.
- 156 8. Woźniak, M.K., Banaszkiwicz, L., Wiergowski, M. *et al.* 2020.
157 Development and validation of a GC–MS/MS method for the determination
158 of 11 amphetamines and 34 synthetic cathinones in whole blood. *Forensic
159 Toxicol*, 38, 42–58.
- 160 9. Wenio I, Bartosiewicz I, Derewiaka D, Dewiszek K, Karniłowicz K. A.
161 2023. Fast Method for Determination of Ethylene Oxide Using Gas
162 Chromatography Coupled with Mass Spectrometry GC-MS/MS. *Applied
163 Sciences*. 13(13):7480. <https://doi.org/10.3390/app13137480>.
- 164 10. Regulation (EC) No 396/2005 of the European Parliament and of the
165 Council of 23 February 2005 on maximum residue levels of pesticides in
166 or on food and feed of plant and animal origin and amending Council
167 Directive 91/414/EEC (OJ L70, 16.3.2005, page 1).